## IMPORTANT NOTE: BIOLOGICAL SPECIMENS FOR METANEPHRINES & CATECHOLS MUST BE KEPT COLD

---

### Background

Catecholamines and other catechols, including DHPG, DOPA and DOPAC, are unstable in most biological fluids at room temperature. The free metanephrines, including normetanephrine and metanephrines, are also unstable. All analytes are subject to auto-oxidation, particularly at ambient temperature and under alkaline conditions. Additionally, catecholamines and catechols in whole blood are susceptible to metabolism by catechol-O-methyltransferase (COMT) in red cells. Thus, while all compounds degrade in whole blood at room temperature, some compounds such as the O-methylated metabolites, can show transient increases during initial periods when COMT might be active. All compounds are somewhat more stable in plasma, but nevertheless will show evidence of degradation when plasma is kept at room temperature. THEREFORE, ALL BLOOD SAMPLES MUST BE PROMPTLY CHILLED AFTER BLOOD COLLECTION AND ALL PLASMA SAMPLES SHOULD BE FROZEN AS SOON AS POSSIBLE AFTER CENTRIFUGATION AND PLASMA IS SEPARATED FROM BLOOD CELLS.

### Methods for Chilling Blood Samples

The best method for chilling blood samples involves their immediate placement on crushed ice. If ice is not available, then other methods can be employed, including the use of cold packs. However, if such methods are employed, care must be taken to ensure that blood is not subject to any freeze effect that may lead to hemolysis. One commercial product suitable for collection of one or two blood samples is detailed below and over the page.

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