

Correct determination of protein concentration in chromaffin cells

Background:

Due to high catecholamine contents in chromaffin cells most protein assays are not suitable for determining the total protein concentration in such cells. One example is given in Fig. 1 where two commercial protein assay reagents, Quant-iT Protein Assay Kit (Invitrogen) and BCA Protein Assay Kit (Pierce), are compared according to their reactivity with norepinephrine. Quant-iT reagent does not react with norepinephrine, whereas with the BCA reagent a concentration-dependent increase in absorbance is evident.




Fig. 1 – Norepinephrine reactivity with Quant-iT Protein Assay reagent (A) and BCA Protein Assay reagent (B)

Material and method:

In order to use Quant-iT Protein Assay Kit (Invitrogen), follow the manufacturer's instructions:

- 1. Equilibrate assay components to room temperature.**
- 2. Make the working solution by diluting Quant-iT™ protein reagent 1:200 in Quant-iT™ protein buffer.**
- 3. Load 200 μ L of working solution in each microplate well.** Diluted Quant-iT™ protein reagent is stable for at least 3 hours at room temperature, protected from light.
- 4. Add 10 μ L of each of the Quant-iT™ protein standards to separate wells and mix well.** Duplicates or triplicates of the standards are recommended.
- 5. Add 1–20 μ L of each unknown protein sample to separate wells and mix well.** Duplicates or triplicates of the unknown samples are recommended.
- 6. Measure fluorescence using microplate reader (excitation/emission maxima \sim 470/570 nm).** The fluorescence signal is stable for 3 hours at room temperature.
- 7. Use a standard curve to determine protein amounts.** For the BSA standards, plot amount vs. fluorescence, and fit a straight line to the data points.